

# GCSE Biology required practical activity: Food tests

## Using qualitative reagents to test for a range of carbohydrates, lipids and proteins

### Risk assessment

- Safety goggles should be worn when carrying out the tests.
- Wash off spills on skin immediately.
- Take care with boiling water.

### Method

#### You are provided with the following:

- food to be tested
- a pestle and mortar
- a stirring rod
- filter funnel and filter paper
- 2 × beaker, 250 ml
- a conical flask
- 2 × test tube
- Benedict's solution
- iodine solution
- kettle for boiling water
- a thermometer
- safety goggles.

#### Read these instructions carefully before you start work.

1. Use a pestle and mortar to grind up a small sample of food.
2. Transfer the ground up food into a small beaker. Then add distilled water.
3. Stir the mixture so that some of the food dissolves in the water.
4. Filter using a funnel with filter paper to obtain as clear a solution as possible.  
The solution should be collected in a conical flask.
5. Half fill a test tube with some of this solution.
6. Add 10 drops of Benedict's solution to the solution in the test tube.
7. Put hot water from a kettle in a beaker. The water should **not** be boiling.  
Put the test tube in the beaker for about five minutes.

8. Note any colour change.

If a reducing sugar (such as glucose) is present, the solution will turn green, yellow, or brick-red. The colour depends on the sugar concentration.

9. Take 5 ml of the solution from the conical flask and put it into a clean test tube.

10. Add a few drops of iodine solution and note any colour change.

If starch is present, you should see a black or blue-black colour appear.

11. Record your results in a table such as the one below.

Name of food tested	Colour produced with Benedict's solution	Colour produced with iodine solution

## 2. Testing for lipids

In this experiment you will test one or more foodstuffs for the presence of lipids (fats).

Learning outcomes
<p>1</p> <p>2</p> <p>3</p> <p>Teachers to add these with particular reference to working scientifically</p>

### Risk assessment:

- Safety goggles should be worn when carrying out the tests.
- Sudan III contains ethanol, which is highly flammable. Keep the solution away from naked flames.
- Wash off spills on skin immediately.

### Method

#### You are provided with the following:

- food to be tested
- a pestle and mortar
- a stirring rod
- 2 × beaker, 250 ml
- a test tube
- Sudan III stain solution.
- safety goggles.

**Read these instructions carefully before you start work.**

1. Use a pestle and mortar to grind up a small sample of food.
2. Transfer the ground up food into a small beaker. Then add distilled water.
3. Stir the mixture so that some of the food dissolves in the water. Do not filter.
4. Half fill a test tube with some of this solution.
5. Add 3 drops of Sudan III stain to the solution in the test tube. Shake gently to mix.
6. If fat is present: a red-stained oil layer will separate out and float on the water surface.

### **3. Testing for proteins**

In this experiment you will test one or more foodstuffs for the presence of protein.

<b>Learning outcomes</b>
<b>1</b>
<b>2</b>
<b>3</b>
<b>Teachers to add these with particular reference to working scientifically</b>

**Risk assessment:**

- Safety goggles should be worn when carrying out the tests.
- Biuret solution contains copper sulphate, which is poisonous, and sodium hydroxide, which is caustic.
- Wash off spills on skin immediately.

**Method**

**You are provided with the following:**

- food to be tested
- a pestle and mortar
- a stirring rod
- a filter funnel and filter paper
- 2 × beaker, 250 ml
- a test tube
- Biuret solution
- safety goggles.

**Read these instructions carefully before you start work.**

1. Use a pestle and mortar to grind up a small sample of food.
2. Transfer the ground up food into a small beaker. Then add distilled water.
3. Stir the mixture so that some of the food dissolves in the water.
4. Filter using a funnel with filter paper to obtain as clear a solution as possible.  
The solution should be collected in a conical flask.
5. Put 2 cm<sup>3</sup> of this solution into a test tube.
6. Add 2 cm<sup>3</sup> of Biuret solution to the solution in the test tube. Shake gently to mix.
7. Note any colour change. Proteins will turn the solution pink or purple.